

Human Crk protein (Recombinant) (His-Tag) (STJP016906)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	25.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1398
Gene Symbol	CRK
UniProt ID	CRK_HUMAN
Immunogen Region	1-204aa
Immunogen Sequence	MAGNFDSEER SSWYWGRLSR QEAVALLQGQ RHGVFLVRDS STSPGDYVLS VSENSRVSHY IINSSGPRPP VPPSPAQPPP GVSPSRLRIG DQEFDSLPL LEFYKIHYLD TTTLIEPVSR SRQSGSVILR QEEAEYVRAL FDFNGNDEED LPFKKGDLR IRDKPEEQWW NAEDSEGKRG MIPVPYVEKY RPASASVSAL IGGR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Crk protein 3µg under reducing conditions stained with Coomassie blue