

Human DNAJB2 protein (Recombinant) (His-Tag) (STJP016928)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.1M NaCl
Molecular Weight	33 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3300
Gene Symbol	DNAJB2
UniProt ID	DNJB2_HUMAN
Immunogen Region	1-277aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMASYYEI LDVPRASAD DIKKAYRRKA LQWHPDKNPD NKEFAEKKFK EVAEAYEVLS DKHKREIYDR YGREGLTGTG TGPSRAEAGS GGPGFTFTFR SPEEVFREFF GSGDPFAELF DDLGPFSELQ NRGSRHSGPF FTFSSSFPGH SDFSSSSF SF SPGAGAFRSV STSTTFVQGR RITTRRIMEN GQERVEVED GQLKSVTI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DNAJB2 protein 3µg under reducing conditions with Coomassie blue stain