

Human 17 beta-HSD14/HSD17B14 protein (Recombinant) (His-Tag) (STJP016937)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	32.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51171
Gene Symbol	HSD17B14
UniProt ID	DHB14_HUMAN
Immunogen Region	1-270aa
Immunogen Sequence	MATG TRYAGKVVVV TGGGRGIGAG IVRAFVNSGA RVVICDKDES GGRALEQELP GAVFILCDVT QEDDVKTLVS ETIRRFGRDL CVVNNAGHHP PPQRPEETSA QGFRQLLELN LLGTYTLTKL ALPYLRKSQG NVINISLVG AIGQAQAVPY VATKGAVTAM TKALALDESP YGVRVNCISP GNIWTP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human 17 beta-HSD14/HSD17B14 protein under reducing conditions