

Human Sirtuin 2/SIRT2 protein (Recombinant) (His-Tag) (STJP016942)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 2mM DTT, 200mM NaCl, 0.5mM EDTA
Molecular Weight	41.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22933
Gene Symbol	SIRT2
UniProt ID	SIR2_HUMAN
Immunogen Region	1-352aa
Immunogen Sequence	MDFLRNLFQ TSLGSKER LDELTLGV ARYMQSERCR RVICLVGAGI STSAGIPDFR SPSTGLYDNL EKYHLPYPEA IFEISYFKKH PEPFFALAKE LYPGQFKPTI CHYFMRLKLD KGLLLRCYTQ NIDTLERIAQ LEQEDLVEAH GTFYTSCHCVS ASCRHEYPLS WMKEKIFSEV TPKCEDCQSL VKPDIVFFGE SLPARF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human SIRT2 protein 3 µg Coomassie blue stain