

Human MP1/LAMTOR3 protein (Recombinant) (His-Tag) (STJP016944)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 5mM DTT, 100mM NaCl
Molecular Weight	15.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8649
Gene Symbol	LAMTOR3
UniProt ID	LTOR3_HUMAN
Immunogen Region	1-124aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MADDLKRFLY KKLPSVEGLH AIVVSDRDGV PVIKVANDNA PEHALRPGFL STFALATDQG SKLGLSKNKS IICYNTYQV VQFNRLPLVV SFIASSANT GLIVSLEKEL APLFEELRQV VEVS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MP1 LAMTOR3 protein under reducing conditions stained with Coomassie blue