

Human HUS1 protein (Recombinant) (His-Tag) (STJP016954)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 40% Glycerol
Molecular Weight	33.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3364
Gene Symbol	HUS1
UniProt ID	HUS1_HUMAN
Immunogen Region	1-280aa
Immunogen Sequence	MKFRAKIVDG ACLNHFTRIS NMIAKLAKTC TLRISPDKLN FILCDKLANG GVSMWCELEQ ENFFNEFQME GVSAENNEIY LETSENLSR ALKTAQNARA LKIKLTNKH F PCLTVSVELL SMSSSSRIVT HDIPIKVIPR KLWKDLQEPV VPDPDVSIL PVLKTMKSVV EKMKNISNHL VIEANLDGEL NLKIETELVC VTTHFK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HUS1 protein under reducing conditions stained with Coomassie blue