

Human MDP1 protein (Recombinant) (His-Tag) (STJP016962)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 100mM NaCl
Molecular Weight	22.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	145553
Gene Symbol	MDP1
UniProt ID	MGDP1_HUMAN
Immunogen Region	1-176aa
Immunogen Sequence	MARLPK LAVFDLDYTL WPFWVDTHVD PPFHKSSDGT VRDRRGQDVR LYPEVPEVLK RLQSLGVPGA AASRTSEIEG ANQLLELFDL FRYFVHREIY PGSKITHFER LQKKTGIPFS QMIFFDDERR NIVDVSKLGV TCIHQNGMN LQTL SQGLET FAKAQTGPLR SSLEESPFEA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MDP1 protein 3µg reducing conditions Coomassie blue stain