

Human RAB6A protein (Recombinant) (His-Tag) (STJP016969)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	26.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5870
Gene Symbol	RAB6A
UniProt ID	RAB6A_HUMAN
Immunogen Region	1-208aa
Immunogen Sequence	MSTGGD FGNPLRKFKL VFLGEQSVGK TSLITRFMYD SFDNTYQATI GIDFLSKTMY LEDRTIRLQL WDTAGQERFR SLIPSYIRDS AAAVVVYDIT NVNSFQQTTK WIDDVRTERG SDVIIMLVGN KTDLADKRQV SIEEGERKAK ELNVMFIETS AKAGYNVKQL FRRVAAALPG MESTQDRSRE DMIDIKLEKP QEQPVS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAB6A protein 3µg under reducing conditions Coomassie blue stain