

Human ROBLD3/LAMTOR2 protein (Recombinant) (His-Tag) (STJP016972)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 10% Glycerol
Molecular Weight	16 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	28956
Gene Symbol	LAMTOR2
UniProt ID	LTOR2_HUMAN
Immunogen Region	1-125aa
Immunogen Sequence	MLRPKA LTQVLSQANT GGVQSTLLN NEGSLAYSG YGDTDARVTA AIASNIWAAY DRNGNQAFNE DNLKFILMDC MEGRVAITRV ANLLLCMYAK ETVGFGMLKA KAQALVQYLE EPLTQVAAS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human ROBLD3 LAMTOR2 protein under reducing conditions with Coomassie blue stain