

Human G3BP1 protein (Recombinant) (His-Tag) (STJP016979)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 100mM NaCl
Molecular Weight	53.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10146
Gene Symbol	G3BP1
UniProt ID	G3BP1_HUMAN
Immunogen Region	1-466aa
Immunogen Sequence	MVMEKPSPLL VGREFVRQYY TLLNQAPDML HRFYGNKSSY VHGGLDNKGK PADAVYGQKE IHRKVMSQNF TNCHTKIRHV DAHATLNDGV VVQVMGLLSN NNQALRRFMQ TFLAPEGVS ANKFYVHNDI FRYQDEVFGG FVTEPQEESE EEEVEPEERQ QTPEVVPDDS GTFYDQAVVS NDMEEHLEEP VAEPEPDPEP EPEQEPVSEI QEEKPEPVLE ETAPEDAQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human G3BP1 protein under reducing conditions stained with Coomassie blue