

Human ARFL2/ARL2 protein (Recombinant) (His-Tag) (STJP016983)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1mM PMSF, 1mM DTT, 40% Glycerol
Molecular Weight	23.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	402
Gene Symbol	ARL2
UniProt ID	ARL2_HUMAN
Immunogen Region	1-184aa
Immunogen Sequence	MGLLTILKKM KQKERELRL MLGLDNAGKT TILKKFNGED IDTISPTLGF NIKTLEHRGF KLNWDVGGQ KSLRSYWRNY FESTDGLIWV VDSADRQRMQ DCQRELQSL VEERLAGATL LIFANKQDLP GALSSNAIRE ALELDSIRSH HWCIGGCSAV TGENLLPGID WLLDDISSRI FTAD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human ARFL2 ARL2 protein under reducing conditions