

Human RAB5B protein (Recombinant) (His-Tag) (STJP016984)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 0.2M NaCl
Molecular Weight	26.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5869
Gene Symbol	RAB5B
UniProt ID	RAB5B_HUMAN
Immunogen Region	1-215aa
Immunogen Sequence	MTSRSTA RPNGQPQASK ICQFKLVLLG ESAVGKSSLV LRFVKGQFHE YQUESTIGAAF LTQSVCLDDT TVKFEIWDTA GQERYHSLAP MYYRGAQAAI VVYDITNQET FARAKTWVKE LQRQASPSIV IALAGNKADL ANKRMVEYEE AQAYADDNSL LFMETSAKTA MNVNDLFLAI AKKLPKSEPP NLGGAAGRSR GVDLHE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAB5B protein 3µg under reducing conditions Coomassie blue stain