

Human Coactosin-like 1/COTL1 protein (Recombinant) (His-Tag) (STJP016996)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	18.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23406
Gene Symbol	COTL1
UniProt ID	COTL1_HUMAN
Immunogen Region	1-142aa
Immunogen Sequence	MATKIDKEAC RAAYNLVRDD GSAVIWVTFK YDGSIVPGE QGAIEYQHFIQ QCTDDVRLFA FVRFTTGDM SKRSKFALIT WIGENVSLQ RAKTGTDKTL VKEVVQNFAK EFVISDRKEL EEDFIKSELK KAGGANYDAQ TE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human COTL1 protein under reducing conditions stained with Coomassie blue