

Human DNAJB6 protein (Recombinant) (His-Tag) (STJP016997)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 30% Glycerol
Molecular Weight	38.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10049
Gene Symbol	DNAJB6
UniProt ID	DNJB6_HUMAN
Immunogen Region	1-326aa
Immunogen Sequence	MVDYYEV LGVQRHASPE DIKKAYRKLA LKWHPDKNPE NKEEAERKFK QVAEAYEVLS DAKKRDIYDK YGKEGLNNGG GGGSHFDSPF EFGFTFRNPD DVFREFFGGR DPFSDFDFED PFEDFFGNRR GPRGSRSRGT GSFFSAFSGF PSFGSGFSSF DTGFTSFGSL GHGGLTSFSS TSFGGSGMGN FKSISTSTKM VNGRKI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human DNAJB6 protein under reducing conditions stained with Coomassie blue