

Human AK4 protein (Recombinant) (His-Tag) (STJP017018)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 2mM DTT, 0.1M NaCl.
Molecular Weight	27.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	205
Gene Symbol	AK4
UniProt ID	KAD4_HUMAN
Immunogen Region	1-223aa
Immunogen Sequence	MASKLLRAVI LGPPGSGKGT VCQRIQNFG LQHLSSGHFL RENIKASTEY GEMAKQYIEK SLLVPDHVIT RLMMSELENR RGQHWLLDGF PRTLQQAEL DKICEVDLVI SLNIPFETLK DRLSRRWIHP PSGRVYNLDF NPPHVGIDD VTGEPLVQQE DDKPEAVAAR LRQYKDVAKP VIELYKSRGV LHQFSGTETN KIWPYV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human AK4 protein under reducing conditions with Coomassie blue stain