

Human NUDT2/APAH1 protein (Recombinant) (His-Tag) (STJP017028)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	19 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	318
Gene Symbol	NUDT2
UniProt ID	AP4A_HUMAN
Immunogen Region	1-147aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MALRACGLII FRRCLIPKVD NNAIEFLLQ ASDGIHHWTP PKGHVEPGED DLETALRETQ EEAGIEAGQL TIEGFKREL NYVARNKPKT VIYWLAEVKD YDVEIRLSHE HQAYRWLGLE EACQLAQFKE MKAALQEGHQ FLCSIEA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human NUDT2
 APAH1 protein under reducing conditions