

Human E2-25K/UBE2K protein (Recombinant) (His-Tag) (STJP017032)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 10% Glycerol, 50mM NaCl
Molecular Weight	26.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3093
Gene Symbol	UBE2K
UniProt ID	UBE2K_HUMAN
Immunogen Region	1-200aa
Immunogen Sequence	MANI AVQRIKREFK EVLKSEETSK NQIKVDLVDE NFTELRGEIA GPPDTPYEGG RYQLEIKIPE TYPFNPPKVR FITKIWHPN I SSVTGAICLD ILKDQWAAAM TLRTVLLSLQ ALLAAAEPPD PQDAVVANQY KQNP EMFKQT ARLWAHVYAG APVSSPEYTK KIENLCAMGF DRNAVIVALS SKSWDV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human E2-25K UBE2K protein under reducing conditions stained with Coomassie blue