

Human eIF-1AY/EIF1AY protein (Recombinant) (His-Tag) (STJP017040)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 100mM NaCl, 20% Glycerol
Molecular Weight	18.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9086
Gene Symbol	EIF1AY
UniProt ID	IF1AY_HUMAN
Immunogen Region	1-144aa
Immunogen Sequence	MPKKNKGK GGKNRRRGKN ENESEKRELV FKEDGQEY AQ VIKMLGNGRL EALCFDGVKR LCHIRGKLRK KVVINTSDII LVGLRDYQDN KADVILKYNA DEARSLKAYG ELPEHAKINE TDTFGPGDDD EIQFDDIGDD DEDIDDI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human eIF-1AY protein 3 µg under reducing conditions
Coomassie blue stain