

Human RAB21 protein (Recombinant) (His-Tag) (STJP017045)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 200mM NaCl, 1mM PMSF
Molecular Weight	27.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23011
Gene Symbol	RAB21
UniProt ID	RAB21_HUMAN
Immunogen Region	18-222aa
Immunogen Sequence	YSF KVVLLGEGCV GKTSLVLRVC ENKFNDKHIT TLQASFLTCK LNIGGKRVNL AIWDTAGQER FHALGPIYYR DSNAILVYD ITDEDSFQKV KNWVKELRKM LGNEICLCIV GNKIDLEKER HVS IQEAESEY AESVGAKHYH TSAKQNKIE ELFLDLCKRM IETAQVDERA KNGSSQPGT ARRGVQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human RAB21 protein under reducing conditions Coomassie blue stain