

Human CIB2 protein (Recombinant) (His-Tag) (STJP017047)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 100mM NaCl, 1mM PMSF
Molecular Weight	23.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10518
Gene Symbol	CIB2
UniProt ID	CIB2_HUMAN
Immunogen Region	1-187aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGNKQTIFTE EQLDNYQDCT FFNKKDILKL HSRFYELAPN LVPMDYRKSP IVHVPMSLII QMPELRENPF KERIVAAFSE DGEGLTFND FVDMFSLCE SAPRELKANY AFKIYDFNTD NFICKEDLEL TLARLTSEL DEEEVVLVCD KVIEEADLDG DGKLGFADE DMIKAPDFL STFHIRI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CIB2 protein under reducing conditions stained with Coomassie blue