

Human IMPA2 protein (Recombinant) (His-Tag) (STJP017051)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT
Molecular Weight	33.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3613
Gene Symbol	IMPA2
UniProt ID	IMPA2_HUMAN
Immunogen Region	1-288aa
Immunogen Sequence	MKPSGEDQAA LAAGPWEECF QAAVQLALRA GQIIRKALTE EKRVTSTKSA ADLVTETDHL VEDLIISELR ERFPSHRFIA EEAAASGAKC VLTHSPTWII DPIDGTCNFV HRFPTVAVSI GFAVRQELEF GVIYHCTEER LYTGRRRGGA FCNGQRLRVS GETDLSKALV LTEIGPKRDP ATLKLFLSNM ERLHAKAHG VRVIGS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human IMPA2 protein stained with Coomassie blue