

Human RPS14 protein (Recombinant) (His-Tag) (STJP017062)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% Glycerol, 100mM NaCl
Molecular Weight	18.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6208
Gene Symbol	RPS14
UniProt ID	RS14_HUMAN
Immunogen Region	1-151aa
Immunogen Sequence	MAPRKGK EKKEEQVISL GPQVAEGENV FGVCHIFASF NDTFVHVTDL SGKETICRVV GGMKVKADRD ESSPYAAMLA AQDVAQRCKE LGITALHIKL RATGGNRTKT PGPGAQSALR ALARSGMKIG RIEDVTPIPS DSTRRKGGRR GRRL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPS14 protein under reducing conditions Coomassie blue stain