

Human NUDT21/CPSF5 protein (Recombinant) (His-Tag) (STJP017072)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 2mM DTT, 20% Glycerol
Molecular Weight	28.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11051
Gene Symbol	NUDT21
UniProt ID	CPSF5_HUMAN
Immunogen Region	1-227aa
Immunogen Sequence	MSVVPNRSQ TGVPRGVTQF GNKYIQQTKP LTLERTINLY PLTNYTFGTK EPLYEKDSSV AARFQRMREE FDKIGMRRTV EGV LIVHEHR LPHVLLQLG TTFKLP GGE LNPGEDEVEG LKRLMTEILG RQDGVLDWV IDDCIGNWWR PNFEPPQYPY IPA HITKPE HKKFLVQLQ EKALFAVPKN YKLVAAPLFE LYDNAP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NUDT21
CPSF5 protein under reducing conditions
Coomassie blue stain