

Human 4E-BP2/EIF4EBP2 protein (Recombinant) (His-Tag) (STJP017101)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	15.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1979
Gene Symbol	EIF4EBP2
UniProt ID	4EBP2_HUMAN
Immunogen Region	1-120aa
Immunogen Sequence	MSSSAGSGHQ PSQSRAIPTR TVAISDAAQL PHDYCTTPGG TLFSTTPGGT RIIYDRKFL DRRNSPMAQT PPCHLPNI PG VTSPGTLIED SKVEVNNLNN LNNHDKRHAV GDDAQFEMDI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human 4E-BP2/EIF4EBP2 protein under reducing conditions