

Human NECAP2 protein (Recombinant) (His-Tag) (STJP017110)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.1M NaCl.
Molecular Weight	30.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55707
Gene Symbol	NECAP2
UniProt ID	NECP2_HUMAN
Immunogen Region	1-263aa
Immunogen Sequence	MEESGYESVL CVKPDVHVYR IPPRATNRYG RAAEWQLDQP SWSGRLRITA KGQMAYIKLE DRTSGELFAQ APVDQFPGTA VESVTSSRY FVIRIEDGNG RRAFIGIGFG DRGDADFNV ALQDHFKWKVQ QCCEFAKQAQ NPDQGPCLDL GFKEGQTIKL NIANMKKKEG AAGNPRVRPA STGGLSLLPP PPGGKTSTLI PPPGEQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human NECAP2 protein SDS-PAGE gel under reducing conditions stained with Coomassie blue