

Human ALKBH3 protein (Recombinant) (His-Tag) (STJP017118)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.1M NaCl.
Molecular Weight	37.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	221120
Gene Symbol	ALKBH3
UniProt ID	ALKB3_HUMAN
Immunogen Region	1-286aa
Immunogen Sequence	M EEKRRRRARVQ GAWAAPVKSQ AIAQPATTAK SHLHQKPGQT WKNKEHHLS D REFVFKEPQQ VVRRAPEPRV IDREGVYEIS LSPTGVSRVC LYPGFVDVKE ADWILEQLCQ DVPWKQRTGI REDITYQQPR LTAWYGELPY TYSRITMEPN PHWHPVLR TL KNRIEENTGH TFNSLLCNLY RNEKDS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human ALKBH3 protein 3µg stained with Coomassie blue