

Human CMPK1 protein (Recombinant) (His-Tag) (STJP017119)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT, 100mM NaCl
Molecular Weight	28.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51727
Gene Symbol	CMPK1
UniProt ID	KCY_HUMAN
Immunogen Region	1-228aa
Immunogen Sequence	MLSRCSRGLL HVLGLSFLQ TRRPILLCSPLMKPLVVFV LGPGGAGKGT QCARIVEKYG YTHLSAGELL RDERKNPDSQ YGELIEKYIK EGKIVPVEIT ISLLKREMDQ TMAANAQKNK FLIDGFPRNQ DNLOGWNKTM DGKADVSVFL FFDCNNEICI ERCLERKSS GRSDDNRESL EKRIQTYLQS TKPIIDLYEE MGKVKK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CMPK1 protein under reducing conditions Coomassie blue stain