

Human Galectin-7/LGALS7 protein (Recombinant) (His-Tag) (STJP017123)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT
Molecular Weight	17.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3963/653499
Gene Symbol	LGALS7.LGALS7B
UniProt ID	LEG7_HUMAN
Immunogen Region	1-136aa
Immunogen Sequence	MSNVPHKSSL PEGIRPGTVL RIRGLVPPNA SRFHVNLLCG EEQGSDAALH FNPRLDTSEV VFNSKEQGSW GREERGPV FQRGQPFV LIIASDDGFKA VVGDAQYHHF RHRLPLARVR LVEVGGDVQL DSVRIF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Galectin-7 LGALS7 protein 3µg reducing conditions Coomassie blue stain