

Human Calcineurin B2/PPP3R2 protein (Recombinant) (His-Tag) (STJP017126)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT
Molecular Weight	22 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5535
Gene Symbol	PPP3R2
UniProt ID	CANB2_HUMAN
Immunogen Region	1-173aa
Immunogen Sequence	MSTMGNEASY PAEMCSHFDN DEIKRLGRRF KKLDLKSGS LSVEEFMSLP ELRHNPLVRR VIDVFDTGDG GEVDFKEFIL GTSQFSVKGD EEQKLRFAS IYDMKDGYI SNGELFQVLK MMVGNNLTDW QLQQLVDKTI IILDKDGDGK ISFEEFSAVV RDLEIHKLV LIV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Calcineurin B2/PPP3R2 protein under reducing conditions