

Human SPOP protein (Recombinant) (His-Tag) (STJP017145)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 50% Glycerol, 0.2M NaCl, 2mM EDTA
Molecular Weight	44.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8405
Gene Symbol	SPOP
UniProt ID	SPOP_HUMAN
Immunogen Region	1-374aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSRVSPPPP AEMSSGPVAE SWCYTQIKVV KFSYMWTTNN FSFCREEMGE VIKSSTFSSG ANDKLKWCLR VNPKGLDEES KDYLSTLYLLL VSCPKSEVRA KFKFSILNAK GEETKAMESQ RAYRFVQGD WGFKKFIRRD FLLDEANGLL PDDKLTLCFE VSVVQDSVNI SGQNTMNMVK VPECRLADEL GGLWENSRT DCCLCVAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SPOP protein under reducing conditions stained with Coomassie blue