

Human SNRPG protein (Recombinant) (His-Tag) (STJP017147)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 50% Glycerol, 0.1M NaCl
Molecular Weight	10.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6637
Gene Symbol	SNRPG
UniProt ID	RUXG_HUMAN
Immunogen Region	1-76aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSAHPPELK KFMDKKLSLK LNGGRHVQGI LRGFDPPFMNLVIDECVEMAT SGQQNNIGMV VIRGNSIIML EALERV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human SNRPG protein under reducing conditions Coomassie stain