

Human ALKBH2 protein (Recombinant) (His-Tag) (STJP017152)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% Glycerol
Molecular Weight	31.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	121642
Gene Symbol	ALKBH2
UniProt ID	ALKB2_HUMAN
Immunogen Region	1-261aa
Immunogen Sequence	MDRFLVKGAQ GLLLRKQEEQ EPTGEEPAVL GGDKESTRKR PRREAPGNNG HSAGPSWRHI RAEGLDCSYT VLFKAEADE IFQELEKEVE YFTGALARVQ VFGKWHVSVPR KQATYGDAGL TYTFSGLTLS PKPWIPVLER IRDHVSGVTG QTFNFVLINR YKDGCDHIGE HRDDERELAP GSPIASVSFG ACRDFVFRHK DSRGKS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human ALKBH2 protein 3µg under reducing conditions
Coomassie blue stain