

Human CCL3/MIP-1 alpha protein (Recombinant) (His-Tag) (STJP017155)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 10mM Sodium Citrate buffer (pH 3.5) containing 20% Glycerol
Molecular Weight	10 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6348
Gene Symbol	CCL3
UniProt ID	CCL3_HUMAN
Immunogen Region	24-92aa
Immunogen Sequence	SLAADTPTA CCFSYTSRQI PQNFIADYFE TSSQCSKPGV IFLTKRSRQV CADPSEEWVQ KYVSDLELSA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human CCL3 MIP-1 alpha protein
 SDS-PAGE under reducing conditions stained
 with Coomassie blue