

Human AIMP3/EEF1E1 protein (Recombinant) (His-Tag) (STJP017158)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	21.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9521
Gene Symbol	EEF1E1
UniProt ID	MCA3_HUMAN
Immunogen Region	1-174aa
Immunogen Sequence	MAAAAEISLL EKSLGLSKGN KYSAQGERQI PVLQTNNGPS LTGLTTIAAH LVKQANKEYL LGSTAEKAI VQQWLEYRVT QVDGHSSKND IHTLLKDLNS YLEDKVYLTG YNFTLADILL YYGLHRFIVD LTVQEKEKYL NVSRLFCHI Q HYPGIRQHLS SVVFIKNRLY TNSH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human AIMP3/EEF1E1 protein under reducing conditions stained with Coomassie blue