

Human NUDT10 protein (Recombinant) (His-Tag) (STJP017165)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	19.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170685
Gene Symbol	NUDT10
UniProt ID	NUDT10_HUMAN
Immunogen Region	1-164aa
Immunogen Sequence	MKCKPNQTRT YDPEGFKKRA ACLCFRSERE DEVLVSSSR YPDRWIVPGG GMEPEEEPPGG AAVREYEEA GVKGKLGRLG GVFEQNQDPK HRTYVYVLTV TELLEDWEDS VSIGRKREWF KVEDAIVLQ CHKPVHAEYL EKLKLGGSPT NGNSMAPSSP DSDP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NUDT10 protein under reducing conditions stained with Coomassie blue