

## Human Tollip protein (Recombinant) (His-Tag) (STJP017202)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.25 mg/mL                                                                                                                                          |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.2M NaCl, 2mM DTT                                                                 |
| Molecular Weight    | 32.4 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                   |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">54472</a>                                                                                                                                                                                                             |
| Gene Symbol        | <a href="#">TOLLIP</a>                                                                                                                                                                                                            |
| UniProt ID         | <a href="#">TOLIP_HUMAN</a>                                                                                                                                                                                                       |
| Immunogen Region   | 1-274aa                                                                                                                                                                                                                           |
| Immunogen Sequence | MATTVSTQQRG PVYIGELPQD FLRITPTQQQ RQVQLDAQAA QQLQYGGAVG TVGRLNITVV QAKLAKNYGM TRMDPYCRLR LGYAVYETPT AHNGAKNPRW NKVIHCTVPP GVDSFYLEIF DERAFSMDDR IAWTHITPE SLRQGKVEDK WYSLSGRQGD DKEGMINLVM SYALLPAAMV MPPQPVLMP TVYQQGVGYV PITGMP |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Tollip protein 3µg under reducing conditions Coomassie blue stain