

Human UCH-L5/UCH37 protein (Recombinant) (His-Tag) (STJP017204)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 30% Glycerol, 200mM NaCl, 0.1mM PMSF, 2mM EDTA
Molecular Weight	39.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51377
Gene Symbol	UCHL5
UniProt ID	UCHL5_HUMAN
Immunogen Region	1-329aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MTGNAGEWCL MESDPGVFTE LIKGFGCRGA QVEEWSLEP ENFEKLKPVH GLIFLKWQP GEEPAGSVVQ DSRLDTIFFA KQVINNACAT QAIVSVLLNC THQDVHLGET LSEFKEFSQS FDAAMKGLAL SNSDVIRQVH NSFARQQMFE FDTKTSAKEE DAFHFVSYVP VNGRLYELDG LREGPIDLGA CNQDDWISAV RPVIEKRI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human UCH-L5/UCH37 protein under reducing conditions