

Human BBS3/ARL6 protein (Recombinant) (His-Tag) (STJP017219)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.2M NaCl, 5mM DTT
Molecular Weight	23.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84100
Gene Symbol	ARL6
UniProt ID	ARL6_HUMAN
Immunogen Region	1-186aa
Immunogen Sequence	MGLLDRLSVL LGLKKKEVHV LCLGLDMSGK TTIINKLKPS NAQSQNILPT IGFSIEKFKS SLSFTVFDM SGQGRYRN LW EHYKKEGQAI IFVIDSSDRL RMVVAKEELD TLLNHPDIKH RRIPILFFAN KMDLRDAVTS VKVSQQLCLE NIKDKPWHIC ASDAIKGEGL QEGVDWLQDQ IQTVKT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human BBS3 ARL6 protein under reducing conditions Coomassie blue stain