

Mouse S100A8 protein (Recombinant) (His-Tag) (STJP017220)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	12.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20201
Gene Symbol	S100a8
UniProt ID	S10A8_MOUSE
Immunogen Region	1-89aa
Immunogen Sequence	MPSELEKALS NLIDVYHNYS NIQGNHHALY KNDFKKMVT ECPQFVQNIN IENLFRELDI NSDNAINFEE FLAMVIKGVV ASHKDSHKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant mouse S100A8 protein under reducing conditions Coomassie blue stain