

## Human AP1S2 protein (Recombinant) (His-Tag) (STJP017221)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.1M NaCl, 2mM DTT
Molecular Weight	20.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">8905</a>
Gene Symbol	<a href="#">AP1S2</a>
UniProt ID	<a href="#">AP1S2_HUMAN</a>
Immunogen Region	1-157aa
Immunogen Sequence	MQFMLLFSRQ GKLRLQKWYV PLSDKEKKKI TRELVQTVLA RPKMKCSFLE WRDLKIVYKR YASLYFCCAI EDQDNELITL EIIHRYVELL DKYFGSVCEL DIIFNFEKAY FILDEFLLGG EVQETSCKNV LKAIEQADLL QEEAETPRSV LEEIGLT

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human AP1S2 protein under reducing conditions Coomassie blue stained