

Human RAB2A protein (Recombinant) (His-Tag) (STJP017233)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 100mM NaCl
Molecular Weight	26.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5862
Gene Symbol	RAB2A
UniProt ID	RAB2A_HUMAN
Immunogen Region	1-212aa
Immunogen Sequence	MAYAYL FKYIIIGDTG VGKSCLLQF TDKRFQPVHD LTIGVEFGAR MITIDGKQIK LQIWDTAGQE SFRSITRSYY RGAAGALLVY DITRRDTFNH LTTWLEDARQ HSNSNMVIML IGNKSDLESR REVKKEEGEA FAREHGLIFM ETSAKTASN EFAFINTAKE IYEKIQEGVF DINNEANGIK IGPQHAATNA THAGNQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human RAB2A protein under reducing conditions Coomassie blue stain