

Human RAB14 protein (Recombinant) (His-Tag) (STJP017242)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	26.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51552
Gene Symbol	RAB14
UniProt ID	RAB14_HUMAN
Immunogen Region	1-215aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMATAPYN YSYIFKYIII GDMGVGKSCSL LHQFTEKKFM ADCPHTIGVE FGTRIIIVSG QKIKLQIWDI AGQERFRAVT RSYRGAAGA LMVYDITRRS TYNHLSSWLT DARNLTNPNT VIILIGNKAD LEAQRDVTYE EAKQFAEENG LLFLEASAKT GENVEDAFLE AAKKIYQNIQ DGSGLDLNAAE SGVQHKPSAP QGGRLTSE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAB14 protein 3µg under reducing conditions