

Human NAT1 protein (Recombinant) (His-Tag) (STJP017245)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.1M NaCl.
Molecular Weight	36.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9
Gene Symbol	NAT1
UniProt ID	ARY1_HUMAN
Immunogen Region	1-290aa
Immunogen Sequence	MDIEAYLERI GYKKS RNKLD LETLTDILQH QIRAVPFENL NIHCGDAMD L GLEAIFDQVV RRNRGGWC LQ VNHLLYWALT TIGFETTMLG GYVYSTPAKK YSTGMIHLL QVTIDGRNYI VDAGFGRSYQ MWQPLELISG KDQPQVPCIF RLTEENGFWY LDQIRREQYI PNEEF LHS DL LEDSKYRKIY SFTLKPRTIE DFESMN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NAT1 protein 3µg under reducing conditions Coomassie blue stain