

e.coli ppa protein (Recombinant) (His-Tag) (STJP017267)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 50mM NaCl
Molecular Weight	21.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-176aa
Immunogen Sequence	MSLLNVPAGK DLPEDIYVVI EIPANADPIK YEIDKESGAL FVDRFMSTAM FYPCNYGYIN HTLSLDGDPV DVLVPTPYPL QPGSVIRCRP VGVLKMTDEA GEDAKLVAVP HSKLSKEYDH IKDVNDLPEL LKAQIAHFFE HYKDLEKGKW VKVEGWENAE AAKAEIVASF ERAKNK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E. coli PPA protein under reducing conditions Coomassie blue stain