

Human KLF4 protein (Recombinant) (His-Cam Tag) (STJP017278)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	58.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9314
Gene Symbol	KLF4
UniProt ID	KLF4_HUMAN
Immunogen Region	11-395aa
Immunogen Sequence	MAHHHHHMA DQLTEEQIAE FKEAFSLFDK DGDGTITTKE LGTVMRS LGQ NPTEAELQDM INEVDADGNG TIDFPEFLTM MARKMKD TDS EEEIREAFRV FDKDGN GYIS AAELRHVMTN LGEKLTDEEV DEMIREADID GDGQVNYEEF VQMMTAKGSM AVSDALLPSF STFASGPAGR EKT LRQAGAP NNRWREELSH MKRLPPVLP G RPYDLAAATV ATDLESGG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human KLF4 protein 3µg under reducing conditions Coomassie blue stain