

Human ARL4/ARL4A protein (Recombinant) (His-Tag) (STJP017283)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	24.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10124
Gene Symbol	ARL4A
UniProt ID	ARL4A_HUMAN
Immunogen Region	1-200aa
Immunogen Sequence	MGNGLSDQTS ILSNLPFSQS FHIVILGLDC AGKTTVLRYL QFNEFVNTVP TKGFNTEKIK VTLGNSKTVT FHFWDVGGQE KLRPLWKSYS RCTDGIVFVV DSVDERMEE AKTELHKITR ISENQGVPL IVANKQDLRN SLSEIEIKL LAMGELSSST PWHLQPTCAI IGDGLKEGLE KLHDMIIKRR KMLRQQKKKR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ARL4 ARL4A protein under reducing conditions with Coomassie blue stain