

Human DCAF7 protein (Recombinant) (His-Tag) (STJP017305)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 1mM DTT
Molecular Weight	33.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10238
Gene Symbol	DCAF7
UniProt ID	DCAF7_HUMAN
Immunogen Region	19-295aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MVIAMNWSVR PDKRFRLALG SFVEEYNNKV QLVGLDEESS EFICRNTFDH PYPTTKLMWI PDTKGVYPDL LATSGDYLRV WRVGETETRL ECLLNKNS DFCAPLTSFD WNEVDPYLLG TSSIDTCTI WGLETGQVLG RVNLVSGHVK TQLIAHDKEV YDIAFSRAGG GRDMFASVGA DGSVRMFDLR HLEHSTIIE DPQHPLL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DCAF7 protein under reducing conditions with Coomassie blue stain