

Human MAF1 protein (Recombinant) (His-Tag) (STJP017324)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	31.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84232
Gene Symbol	MAF1
UniProt ID	MAF1_HUMAN
Immunogen Region	1-256aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MKLLENSFE AINSQTVET GDAHIIGRIE SYCKMAGDD KHMFKQFCQE GQPHVLEALS PPQTSGLSPS RLSKSQGGEE EGPLSDKCSR KTLFYLIATL NESFRPDYDF STARSHEFSR EPSLSWVNA VNCSLFSAVR EDFKDLKQL WNAVDEEICL AECDIYSYNP DLSDPFGED GSLWSFNFF YNKRLKRIVF FSCRSISG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human MAF1 protein under reducing conditions stained with Coomassie blue