

Human HSPB3 protein (Recombinant) (His-Tag) (STJP017329)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) 5% Glycerol, 1M urea.
Molecular Weight	19.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8988
Gene Symbol	HSPB3
UniProt ID	HSPB3_HUMAN
Immunogen Region	1-150aa
Immunogen Sequence	MGSSHHHHHH SGLVPRGSH MAKIILRHLI EIPVRYQEEF EARGLEDCLRL DHALYALPGP TIVDLRKTRA AQSPPVDSAA ETPPREGKSH FQILLDVVQF LPEDIIIQTF EGWLLIKAQH GTRMDEHGFI SRSFTRQYKL PDGVEIKDLS AVLCHDGILV VEVKDPVGTK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human HSPB3 protein under reducing conditions with Coomassie blue stain