

Human ARF3 protein (Recombinant) (His-Tag) (STJP017342)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT
Molecular Weight	22.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	377
Gene Symbol	ARF3
UniProt ID	ARF3_HUMAN
Immunogen Region	1-181aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGNIFGNLLK SLIGKKEMRI LMVGLDAAGK TTILYKCLKG EIVTTIPTIG FNVETVEYKN ISFTVWDVGG QDKIRPLWRH YFQNTQGLIF VVDSNDRERV NEAREELMRM LAEDEL RDAV LLVFANKQDL PNAMNAEIT DKLGHLRLH RNWYIQATCA TSGDGLYEGL DWLANQLKNK K

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ARF3 protein 3µg under reducing conditions Coomassie stain